



Technical Regulatory Topic

New ASTM Standard Practice for Extraction of Particulate Matter from the Single-Use Systems

Document Number: USTR3430
Revision Number: 1.0
Date: 07-DEC-2020
Author: Patrick Evrard

Introduction

In May 2020, ASTM[♦] International published a new standard practice on a method to extract particles from single-use systems (SUS). This document, “E3230-20: Standard Practice for Extraction of Particulate Matter from the Surfaces of Single-Use Components and Assemblies Designed for Use in Biopharmaceutical Manufacturing”¹, is the outcome of several years of collaboration inside the ASTM Committee E55 on Manufacture of Pharmaceutical and Biopharmaceutical Products.

Scope of the ASTM Standard Practice

This standard practice provides useful recommendations on validating the effectiveness of an extraction method used to collect loose particles from the surface of single-use components and assemblies. The approach was inspired by an automotive industry standard² used for parts which require critical cleanliness. The guidance on the method for counting the particles is not provided.

1 Can be bought on ASTM International website:
<https://www.astm.org/Standards/E3230.htm> (accessed December 7, 2020).

2 ISO 16232:2018

Content of the ASTM Standard Practice

The ASTM practice contains useful guidance on the following:

- types of extraction method;
- the key elements to consider when defining a liquid extraction method to collect particles from SUS;
- different parameters that can potentially impact the effectiveness of the extraction method;
- tests required to verify that the extraction is done in adequate conditions, for example a background control test;
- qualifying the effectiveness of the liquid extraction method;
- acceptance criteria based on a declining criterion, in which the extraction procedure must achieve a minimum 90% particle removal effectiveness during qualification;
- troubleshooting and optimization of the extraction method during the development and qualification steps.

This standard provides excellent guidance to verify that the extraction method applied is actually effective.



Corporate Headquarters

Port Washington, NY, USA
+1.800.717.7255 toll free (USA)
+1.516.484.5400 phone

European Headquarters

Fribourg, Switzerland
+41 (0) 26 350 53 00 phone

Asia-Pacific Headquarters

Singapore
+65 6389 6500 phone

Visit us on the Web at www.pall.com/biotech
Contact us at www.pall.com/contact

International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States, and Venezuela. Distributors in all major industrial areas of the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2020 Pall Corporation. The Pall logo and Pall are trademarks of Pall Corporation.
® indicates a trademark registered in the USA and TM indicates a common law trademark.
Filtration. Separation. Solution is a service mark of Pall Corporation.
♦ ASTM is a trademark of American Society Corporation for Testing and Materials.