

Subject: Notification of Pall Corporation's Allegro™ single-use systems with component article containing an SVHC (Substance of Very High Concern) in excess of 0.1%

Reference: European REACH Regulations (EC 1907/2006), plus UK implementation statues, and amendments in respect of ECHA SVHC levels in articles and component articles supplied by Pall Corporation.

Silicone based materials are heavily used throughout the biotech industry as elastomeric seals, gaskets, and tubing. On June 27, 2018 decision ED/61/2018 was published by the European Chemicals Agency (ECHA) which added 3 cyclo-siloxane compounds to the European REACH Annex XIV 'Candidate list of substances of very high concern'. These cyclo-siloxane compounds, D4, D5 and D6 being commonly used in the production of silicone-based materials.

Name	CAS Number
Octamethylcyclotetrasiloxane (D4)	556-67-2
Decamethylcyclopentasiloxane (D5)	541-02-6
Dodecamethylcyclohexasiloxane (D6)	540-97-6

Allegro single-use systems are considered 'complex objects', as provided to our customers in the European Economic Area (EEA) and United Kingdom (UK). Please be aware that some Allegro single-use systems incorporate 'component articles' where our suppliers have advised the component contains a silicone material with one of more of D4, D5, D6 siloxanes above the 0.1% level. The identified component articles are certain brands of silicone tubing, silicone O-rings, silicone gaskets, connectors, and aseptic connectors (all non-Pall manufactured) with silicone gaskets and silicone O-rings.

As of January 5, 2021 Pall Corporation has provided entries for those Allegro single-use systems employing component articles with a REACH SVHC above 0.1% in the ECHA SCIP database – see <https://echa.europa.eu/scip>. This is in line with the requirements of Waste Framework Directive 2008/98/EC Article 9(1)(i) as amended by 2018/851/EU.

Pall Corporation has also notified Allegro customers of this action in a customer notification letter, reference number CN-RPII-BURQKN dated November 16, 2020. The customer notification letter also advised that Allegro single-use system labels would be updated to include a link to a regulatory letter (this letter). The label update action has also been completed.

For your reference; Addendum 1 to this letter contains a listing of the ECHA SCIP database entry numbers made by Pall Corporation alongside a description of the applicable component articles in the "Primary Article Identifier" field of SCIP. To check if your specific Allegro single-use system part number is covered by an entry, please consult the field marked "Other article Identifier" within the ECHA SCIP database, where the impacted Allegro single-use system part numbers (subject to that entry) have been listed. Alternatively, you can consult your Allegro single-use system drawing(s) for the component part numbers used in that system(s).

Addendum 2 to this letter contains a listing of all component part numbers which have been entered into the ECHA SCIP database. Please note that Addendum 2 is complete as of the date of this letter but is subject to continuous review as new supplier information is received or new components require additional entries into the ECHA SCIP database. We will update this letter on a periodic basis, but you should consider the ECHA SCIP database as the primary and most up-to-date source of information.



Pall Corporation wishes to advise you that the current REACH regulation does not impose restrictions on the use of silicone-based articles, containing D4, D5, D6 siloxanes, in biopharmaceutical production use. However, the REACH regulation does obligate all manufactures or importers of such articles, operating within the EEA or UK, to communicate the presence of these substances if above the 0.1% level to their customer. Pall Corporation has taken steps to notify with respect to Allegro single-use systems, in order to comply with this aspect of the REACH regulation.

Prepared by: Pall Quality Assurance and Regulatory Affairs

Date of Issue: November 3, 2021

To the best of our knowledge this information is accurate as of the date of issuance. However, these statements are subject to change as new information becomes available. We recommend that you periodically confirm this information.

© Copyright 2021 Pall Corporation. All rights reserved.

**Addendum 1:
Pall Corporation SCIP entry tracker**

SCIP Number	Description in SCIP
d7f85d9a-987b-4d3f-8ad2-8322d242fb8e	Pall Customised Disposable System with Silicone Tube Type 4 (Combination 1)
fa4d1c88-e9ca-4da7-9ade-864aa027549e	Pall Customised Disposable System with Silicone Tube Type 4 and O Ring (Combination 1, 13)
d7c9f2cd-8256-428a-a9fb-755249f9cbd3	Pall Customised Disposable System with Silicone Tube Type 4, O Ring and Silicone Tube Type 3 (Combination 1, 13, 14)
70c3c338-a14c-4534-a3f0-47f071146bde	Pall Customised Disposable System with Silicone Tube Type 4 and Silicone Tube Type 3 (Combination 1, 14)
9b8fa841-3fe3-469c-abb5-256cc1b8a52b	Pall Customised Disposable System with Silicone Tube Type 4 and Silicone Tube Type 3 and Silicone Tube Type 5 (Combination 1, 14, 22)
ab35879c-cba5-4f9d-93ee-0b3e66e80244	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Assembly (Combination 1, 16)
8460df40-b460-4ffc-a2a5-977192a45ac8	Pall Customised Disposable System with Silicone Tube Type 4, Stopper assembly (Combination 1, 18)
69cef7f0-21e5-44f6-9ad0-2f86f961e200	Pall Customised Disposable System with Silicone Tube Type 4, Check Valve (Combination 1, 21)
56f43688-88d9-4d8b-b30d-6aeaeab415ce	Pall Customised Disposable System with Silicone Tube Type 4 and Silicone Tube Type 5 (Combination 1, 22)
05fff30b-93bc-4464-b6cd-1bb36daa0a9f	Pall Customised Disposable System with Silicone Tube Type 4 and Sanitary Connector (Combination 1, 4)
f6a95e4a-1609-4b8b-9438-849d68841e75	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector and O-ring (Combination 1, 4, 13)
ae60916c-f6d0-4b8c-8707-de5cbbbab02f	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector, O ring and Silicone Tube Type 3 (combination 1, 4, 13, 14)
d3ed7223-7f57-4fed-8370-25395f529809	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector and Silicone Tube Type 3 (Combination 1, 4, 14)
b32e676d-bee2-4fd7-9d03-8817c2108236	Pall Customised Disposable System with Silicone Tube Type 4 and Sanitary Connector and Silicone Tube Type 5 (Combination 1, 4, 22)
40eaf4f9-5c97-4aad-a83e-cff98faa12ba	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector and Silicone Tube Type 2 (Combination 1, 4, 8)
6b19c9e7-fb97-4c1a-bb63-a25bcb5d9281	Pall Customised Disposable System with Silicone Tube type 4, Sanitary Connector, Silicone Tube Type 2 and O ring (Combination 1, 4, 8, 13)
768b1467-284f-45ed-9b50-4fc21e5aa551	Pall Customised Disposable System with Silicone Tube type 4, Sanitary Connector, Silicone Tube Type 2, O ring and Silicone Tube Type 3 (Combination 1, 4, 8, 13, 14)
5c1631b5-8b2a-4984-8b2d-c5883d145096	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector and Silicone Tube Type 2 and Silicone Tube Type 5 (Combination 1, 4, 8, 22)

9e729fc0-9c59-4ec8-9123-c45d3f7c0657	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector, Silicone Tube Type 2 and Braided Tubing (Combination 1, 4, 8, 9)
4683033f-ad18-46ac-a1b1-4fd94c839c26	Pall Customised Disposable System with Silicone Tube type 4, Sanitary Connector, Silicone Tube Type 2, Braided Tubing and O ring (Combination 1, 4, 8, 9, 13)
1d075daf-ee1f-44bb-9989-536bf911079a	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector, Silicone Tube Type 2 and Braided Tubing and Silicone Tube Type 5 (Combination 1, 4, 8, 9, 22)
c8d92b3c-82f5-4cd9-bfa6-87f548819772	Pall Customised Disposable System with Silicone Tube Type 4, Sanitary Connector and Braided Tubing (Combination 1, 4, 9)
19dc1d9d-b47d-460b-9016-1a792d2395ab	Pall Customised Disposable System with Silicone Tube Type 4 and Silicone Tube Type 2 (Combination 1, 8)
98d9b314-ab9a-4db4-8892-277d340aeef	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Type 2 and O Ring (Combination 1, 8, 13)
ae8e946e-50e8-4fbf-a65e-17e5f19c02a6	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Type 2 and O Ring and Silicone Tube Type 5 (Combination 1, 8, 13, 22)
2c9c835d-ce62-42ed-a93a-68c8efd80b34	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Type 2 and Silicone Tube Type 3 (Combination 1, 8, 14)
89e2e112-f092-401c-a45b-bbee218eaf65	Pall Customised Disposable System with Silicone Tube Type 4 and Silicone Tube Type 2 and Silicone Tube Type 5(Combination 1, 8, 22)
b4894862-bc2f-48fa-9da3-b689b9e18423	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Type 2 and Braided Tubing (Combination 1, 8, 9)
8a70b5d8-66b8-4aec-8e31-d3bed0b822bc	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Type 2, Braided Tubing and Silicone Tube Type 3 (Combination 1, 8, 9, 14)
26c0ef72-4ea2-465b-8bf9-230a93460e05	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Type 2, Braided Tubing and Silicone Tube Type 3 and Silicone Tube Type 5(Combination 1, 8, 9, 14, 22)
d556f946-af0a-4bd5-b125-f4375e8afc9d	Pall Customised Disposable System with Silicone Tube Type 4, Silicone Tube Type 2 and Braided Tubing and Silicone Tube Type 5 (Combination 1, 8, 9, 22)
47ff586e-3836-4808-9604-c102ba76b303	Pall Customised Disposable System with Silicone Tube Type 4 and Braided Tubing (Combination 1, 9)
b7559351-5de9-48f1-926a-b7d11d775c57	Pall Customised Disposable System with Silicone Tube Type 4 and Braided Tubing With O ring and Silicone Tube Type 3(Combination 1, 9, 13, 14)
951f72a2-e930-42f9-928b-143c14329b4d	Pall Customised Disposable System with Tube Assembly (Combination 12)
dc6ab6ce-0b5c-44d9-9526-2ca63d596353	Pall Customised Disposable System with O Ring (Combination 13)
00043bc8-342c-474d-bf2b-bd13eaedc982	Pall Customised Disposable System with O ring and Silicone Tube Type 3 (Combination 13, 14)



259cd57a-513b-4cd9-9e82-d43cf3bb74c3	Pall Customised Disposable System with Silicone Tube Type 3 (Combination 14)
f40408d8-41c4-41f2-bfa0-d20ce9021d32	Pall Customised Disposable System with Liner (Combination 15)
6a584396-421a-4bf7-850e-b061b12420fc	Pall Customised Disposable System with Silicone Tube Assembly (Combination 16)
81bb93e8-9307-44e0-ae0b-a9e7548acdb2	Pall Customised Disposable System with Caps (Combination 17)
ee4505df-174a-4a0e-8195-9c18789bee3d	Pall Customised Disposable System with Silicone Gasket (Combination 19)
fa10832f-d951-4f88-934c-7949cf6dcd72	Pall Customised Disposable System with Aseptic, Quick Connector (Combination 2)
a275a806-f16c-4e77-ae84-34c39f51bee1	Pall Customised Disposable System with Aseptic, Quick Connector and Silicone Tube Type 2 (Combination 2, 8)
c8df70a7-8567-4daa-bc31-600c2c94e6b1	Pall Customised Disposable System with Aseptic, Quick Connector, Silicone Tube Type 2 and O ring (Combination 2, 8, 13)
bafd02d8-bc4d-45ec-b6c0-249e2c4fee7c	Pall Customised Disposable System with Aseptic, Quick Connector and Silicone Tube Type 2 and Braided Tubing (Combination 2, 8, 9)
e602f629-dde4-49a4-b192-d1cc1f27ace1	Pall Customised Disposable System with Aseptic, Quick Connector, Silicone Tube Type 2, Braided Tubing and O ring (Combination 2, 8, 9, 13)
d5ff39c5-0e11-4999-80a6-1213d55726d6	Pall Customised Disposable System with Aseptic, Quick Connector and Braided Tubing (Combination 2, 9)
dd672679-ad5e-466f-ab55-abfa8931afd2	Pall Customised Disposable System with Aseptic, Quick Connector, Braided Tubing and O ring (Combination 2, 9, 13)
376a2d02-88fa-4727-8806-3204f148e8e8	Pall Customised Disposable System with Silicone Gasket Type 2 (Combination 20)
8914c7d1-8f43-4728-a507-f0a235fc9093	Pall Customised Disposable System with Silicone Check Valve (Combination 21)
5fe185f2-8d63-47cd-8840-68edbf1bf1a5	Pall Customised Disposable System with Silicone Tube Type 5 (Combination 22)
01d52c5c-84e8-4c22-ae50-1713d8a23faa	Pall Customised Disposable System with Sanitary Connector (Combination 4)
49a68f1e-2081-4299-b980-50ae90ab45fe	Pall Customised Disposable System with Sanitary Connector and O ring (Combination 4, 13)
3e0b70c1-5cdf-4431-a073-24155b0e2a35	Pall Customised Disposable System with Sanitary Connector and Silicone Tube Type 3 (Combination 4, 14)
e37811d4-db15-465d-b383-f424b5412f6a	Pall Customised Disposable System with Sanitary Connector and Silicone Tube Type 3 and Silicone Tube Type 5 (Combination 4, 14, 22)
a6e8f351-58f4-4670-9225-c103a21342c6	Pall Customised Disposable System with Sanitary Connector and Liner (Combination 4, 15)
d619d71c-4288-4070-a93b-bf113bdb7fdf	Pall Customised Disposable System with Sanitary Connector and Silicone Tube Type 5 (Combination 4, 22)
dc8d3034-afb0-452c-b054-dc9fda07b531	Pall Customised Disposable System with Sanitary Connector and Silicone Tube Type 2 (Combination 4, 8)



1527e423-c5bc-4e40-b4a0-23a918125590	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2 and O ring (Combination 4, 8, 13)
d5f008a4-647e-46ba-a6ee-5d956452d4a2	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2, O ring and Silicone Tube Type 3 (Combination 4, 8, 13, 14)
1c78245d-9135-4702-b2c8-e4b28350b38c	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2 and Silicone Tube Type 3 (Combination 4, 8, 14)
0c2cbe17-cedd-4caa-a335-3ce5949076a3	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2 and Liner (Combination 4, 8, 15)
67c537d3-2eaa-41c3-b74c-5c85b2b54aae	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2 and Caps (Combination 4, 8, 17)
9e7bfbeb-407e-49aa-8ac0-eea932711701	Pall Customised Disposable System with Sanitary Connector and Silicone Tube Type 2 and Silicone Tube Type 5 (Combination 4, 8, 22)
c5a5b541-a136-4004-bafc-992242a26a9d	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2 and Braided Tubing (Combination 4, 8, 9)
61785cd7-24c3-4c92-92b3-e44335178e20	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2, Braided Tubing and O ring (Combination 4, 8, 9, 13)
0c8512ae-786a-48cd-8192-7124e951848a	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2, Braided Tubing and Silicone Tube Type 3 (combination 4, 8, 9, 14)
bec39b25-1a3a-46b3-8070-931464bb4a44	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2, Braided Tubing and Liner (Combination 4, 8, 9, 15)
7d16be1e-985f-44c4-b4df-26ea74d57e56	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2, Braided tubing and Caps (Combination 4, 8, 9, 17)
199c9634-8b83-4d64-99f2-79538214ae90	Pall Customised Disposable System with Sanitary Connector, Silicone Tube Type 2, Braided tubing, and Silicone Gasket Type 2 (Combination 4, 8, 9, 20)
bc4ee01f-30b9-4c69-b6eb-b70322aab807	Pall Customised Disposable System with Sanitary Connector and Braided Tubing (Combination 4, 9)
0566886b-0e73-4ec6-81ba-e2ea8ce103e7	Pall Customised Disposable System with Sanitary Connector, Braided Tubing and O ring (Combination 4, 9, 13)
23bef0ec-7848-4141-b20a-17bc5c455849	Pall Customised Disposable System with Sanitary Connector, Braided Tubing and Liner (Combination 4, 9, 15)
dc7291e7-6312-4632-afef-df5c43944ee5	Pall Customised Disposable System with Silicone Tube Type 2 (Combination 8)
e5b49a47-51c5-4109-abbce3e922c29153	Pall Customised Disposable System with Silicone Tube Type 2 and Tube Type 1 (Combination 8, 10)
866cbcb6-999d-4767-b4ffe0c8e5b2f90	Pall Customised Disposable System with Silicone Tube Type 2 and Tube Assembly (Combination 8, 12)
338aea43-02df-418b-9c55-fe3837d96f26	Pall Customised Disposable System with Silicone Tube Type 2 and O ring (Combination 8, 13)
3b5d3f6f-c81b-467a-afed-2901a472ccca	Pall Customised Disposable System with Silicone Tube Type 2, O ring and Silicone Tube Type 3 (Combination 8, 13,14)
0ef1fc78-0847-452e-ba00-6c9a8ed8835f	Pall Customised Disposable System with Silicone Tube Type 2, O ring and Silicone Gasket Type 2 (Combination 8, 13,20)



676e21df-20f8-4c61-8b03-1c3a2fcf5297	Pall Customised Disposable System with Silicone Tube Type 2 and Silicone Tube Type 3 (Combination 8, 14)
e94674a3-5e3c-43e8-baa5-18825242cd4a	Pall Customised Disposable System with Silicone Tube Type 2 and Liner (Combination 8, 15)
1c3e17f5-a6f9-4a56-95e6-ee2785e9efad	Pall Customised Disposable System with Silicone Tube Type 2 and Silicone Tube Assembly (Combination 8, 16)
8afce207-118f-4fde-a26c-7eb0a1340f57	Pall Customised Disposable System with Silicone Tube Type 2 and Caps (Combination 8, 17)
c087d0fa-b2a8-4201-b390-6ef873f5449f	Pall Customised Disposable System with Silicone Tube Type 2 and Stopper Assembly (Combination 8, 18)
6e53e609-bc4c-4f74-a0a7-e35f0225013c	Pall Customised Disposable System with Silicone Tube Type 2, Stopper assembly, and Caps (Combination 8, 18, 17)
1ce291f3-f1de-421a-87c4-9fb9cbefcfff	Pall Customised Disposable System with Silicone Tube Type 2 and Silicone Gasket (Combination 8, 19)
a3edeb5e-7009-46eb-9a34-f0e2396d8245	Pall Customised Disposable System with Silicone Tube Type 2, and Silicone Gasket Type 2 (Combination 8, 20)
b8d5736c-8fdc-40a6-acef-b39a761b0aac	Pall Customised Disposable System with Silicone Tube Type 2, and Check Valve (Combination 8,21)
c8c7ab27-6323-42b5-8060-5e6608eca238	Pall Customised Disposable System with Silicone Tube Type 2 and Silicone Tube Type 5 (Combination 8, 22)
cd83fe23-4bbf-4c47-9066-8b2053c9f159	Pall Customised Disposable System with Silicone Tube Type 2 and Braided Tubing (Combination 8, 9)
bdbf1a43-650b-4204-a201-cabc00040648	Pall Customised Disposable System with Silicone Tube Type 2, Braided Tubing and Tube Assembly (Combination 8, 9, 12)
d19075ae-7e70-464c-868f-61b37cfea0b0	Pall Customised Disposable System with Silicone Tube Type 2, Braided Tubing and O ring (combination 8, 9, 13)
17404956-2d84-4871-ba39-dd035220de87	Pall Customised Disposable System with Silicone Tube Type 2, Braided Tubing and Silicone Tube Type 3 (Combination 8, 9, 14)
443a6bd7-89b5-4b25-8670-a0f5fb577615	Pall Customised Disposable System with Silicone Tube Type 2, Braided Tubing and Liner (Combination 8, 9, 15)
de786e2b-06d1-4ec9-a4b3-a821e2f94b1b	Pall Customised Disposable System with Silicone Tube Type 2, Braided tubing and Caps (Combination 8, 9, 17)
ecb4866e-c2af-4628-8bb5-29d3caaa4628	Pall Customised Disposable System with Silicone Tube Type 2, Braided Tubing and Silicone Gasket (Combination 8, 9, 19)
3fbea93b-62dd-4bd8-b37d-b74f0f4a2968	Pall Customised Disposable System with Silicone Tube Type 2, Braided Tubing and Silicone Gasket type 2 (Combination 8, 9, 20)
65f44543-7903-42d7-96df-58d67c7213da	Pall Customised Disposable System with Silicone Tube Type 2 and Braided Tubing and Silicone Tube Type 5 (Combination 8, 9, 22)
872ba668-fa94-46cf-9640-cbbbee53ab37	Pall Customised Disposable System with Braided Tubing (Combination 9)
92fdea5a-04ed-4b94-a2da-3a521b3d7b2d	Pall Customised Disposable System with Braided Tubing and O ring (Combination 9, 13)
46b5f6d5-e3cc-4647-9589-7ae3f05e7bfb	Pall Customised Disposable System with Braided Tubing and Silicone Tube Type 3 (Combination 9, 14)



PALL CORPORATION

aac8f0e1-0d73-456f-abb3-e7dad07d4097	Pall Customised Disposable System with Braided Tubing and Liner (Combination 9, 15)
aec606a-664c-4c18-8ccb-8d492773f9eb	Pall Customised Disposable System with Braided tubing and Caps (Combination 9, 17)
c1c583d1-07b8-4275-ad1a-d83150a197ee	Pall Customised Disposable System with Braided Tubing and Silicone Gasket type 2 (Combination 9, 20)

**Addendum 2:
Component part numbers listed in ECHA SCIP “Other article Identifier” (identified as containing Cyclosiloxane compounds above 0.1%)**

2100554M	250664000	FY02910	SC00982	FY03041	SC01391	SC01813	SC01166	SC01519
NL03203	250665000	NL03496	SC00988	FY03042	SC01753	SC01113	SC01168	SC01917
NL03204	FY02034	SC00113	SC01011	851244000	SC01752	SC01754	SC00357	FY02589
NL03205	FY02035	SC00114	SC01597	851248000	SC01684	SC01751	SC00935	SC02271
NL03206	FY02036	SC00802	SC01519	NL00259	SC01431	SC01344	SC02127	SC01933
NL03207	FY02037	SC01003	NL04218	NL03878	SC01426	SC01345	FY02128	2100618M
NL03415	FY02038	SC01473	NL04223	NL03879	SC01371	SC01346	FY02160	2100619M
NL03498	FY02048	SC01474	NL04229	NL03959	SC01130	SC01347	FY02883	2100620M
NL03508	FY02049	NL04220	NL05187	NL04536	SC01112	SC01348	SC00447	2100621M
NL03833	FY02058	NL04221	SC00811	NL04553	SC01127	SC01734	SC00596	2100622M
NL04064	FY02065	NL04222	SC01096	NL04699	SC00822	SC01735	STC2020900	2100630M
NL04770	FY02089	NL04224	SC01153	NL04701	SC01388	SC01736	SC00236	1101908
NL04777	FY02192	NL04228	SC01538	NL04746	SC00818	SC01737	SC01048	1102399
NL04817	FY02193	NL08214	SC01745	NL05125	SC01111	SC01738	SC01122	1102638
NL05191	FY02202	SC00035	SC01537	NL05178	SC01102	SC02070	851245000	1104545
SC00219	FY02203	SC00036	SC01539	SC00187	SC00820	SC02071	51606	705331
SC00297	FY02204	SC00037	SC01540	SC00826	SC00819	SC02072	SC01155	701563
SC00351	FY02205	SC00038	SC01541	SC01390	SC00821	SC02073	SC01156	SC01712
SC00964	FY02425	SC00039	SC01624	SC01375	SC01877	SC02130	SC01261	SC01716
701560	FY02594	SC00475	SC01663	SC01374	SC00823	SC02120	SC01262	SC01725
701561	FY02622	SC00836	SC01664	SC01373	SC00825	SC02155	SC01281	SC01726
701562	FY02697	SC00957	SC01814	SC01389	SC00827	SC02170	SC01284	SC02140
SC01815	SC00981	SC00974	SC01126	SC00824	SC01017	SC01883	SC01285	SC01881
SC01559	SC02340	SC01744	SC01835					